
FX 100 Fixed Line Frequency



100 mA Fixed Diagnostic X Ray Machine (LINE FREQUENCY)







Tubehead



Controller





TECHNICAL SPECIFICATION

1.Radiography:100 mA at 100 Kvp

2.Tube current : F(0-3 mA), 25 mA, 50 mA, 75 mA & 100 mA

3.Tube voltage : 45-100 kvp (15 Steps)

4.Timer : 0.02-5 Seconds (24 Steps

5.X-ray tube: Indian / Imported Stationery Anod 6.Collimator: Manual Light Beam Diaphragm

7.Power supply :Single Phase 230 + 10 V,50 Hz., Line Resistance 0.4 ohm Maximum

8. Rectification: Full Wave Rectified Using Silicon Rectifiers



- Rails for horizontal movement of the column
- Tube stand and tube head rotatable in 360°
- User friendly design
- Compatible with horizontal as well as multi cone table



AERB APPROVED



X Ray India Private Limited

Contact Us: +91-141-2983862, +91-7073577706 Email: info@xrayindia.net Web: www.xrayindia.net









